

SmartPeg Reference List

valid from November 2014

IMPLANT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
3i	Osseotite/Nanotite External Hex	3.75, 4.0, 5.0, 6.0	Type 1	100353
	Osseotite/Nanotite External Hex	XP 3/4, 4/5, 5/6	Type 1	100353
	Osseotite/Nanotite External Hex MicroMiniplant	3.25	Type 1	100353
	Osseotite/Nanotite Certain MicroMiniplant	3.25	Type 45	100472
	Osseotite/Nanotite Certain Prevail	4/3	Type 45	100472
	Osseotite/Nanotite Certain	4.0, 5.0, 6.0	Type 15	100386
	Osseotite/Nanotite Certain Prevail	5/4, 4/5/4, 5/6/5	Type 15	100386
	Osseotite Certain	XP 4/5, 5/6	Type 15	100386
	Osseotite Tapered Certain	3.25	Type 45	100472
	Osseotite Tapered Certain	4.0	Type 15	100386
	T3 Implant Internal Connection		Type 15	100386
	T3 Implant External Connection		Type 1	100353
	A.B. DENTAL DEVICES	Screw Type, Conical Groovy	3.75	Type 32
4.5			Type 27	100431
ADDENTA	Bifasico	4.0	Type 14	100374
ADIN DENTAL IMPLANT SYSTEMS	One	3.3, 4.25	Type 17	100390
	Swell	3.3 - 6.2	Type 49	100480
	Touareg, Touareg-S, Touareg-X		Type 49	100480
	Touareg CloseFit	RP 3.5	Type 60	100543
	Touareg X-1330	3.0	Type 62	100545
ALPHA BIO IMPLANT	ATI, SPI		Type 26	100425
	ATID	3.75	Type 26	100425
ALFA-GATE	Bioactive	3.75	Type 35	100446
ALPHATECH		3.4 - 7.0	Type 46	100477
ANKYLOS	C/X	A, B, C	Type 16	100388
ANTHOGRYR	Anthofit OI Tapered, OI Straight		Type 18	100398
	Anthofit HE Straight		Type 1	100353
	Axiom	3.4-5.2	Type 38	100455
	Ossfit		Type 4	100350
ARDS IMPLANTS	Smart, Premium, Classic		Type 27	100431
ASTRA TECH	MicroThread	3.5	Type 5	100376
		4.0	Type 6	100378
		4.5, 5.0	Type 7	100380
	Osseospeed Aqua	3.5 S, 4.0 S	Type 38	100455
	OsseoSpeed EV	3.0S Green	Type 50	100481
		3.6S Purple	Type 38	100455
		4.2S Yellow	Type 49	100480
		4.8S Blue	Type 47	100478
		5.4S Brown	Type 40	100462
	Osseospeed Lilac	4.5, 5.0, 5.0 S	Type 7	100380
	Osseospeed TX Aqua	3.5 S, 4.0 S	Type 38	100455
	Osseospeed TX Lilac	4.5, 5.0, 5.0 S	Type 7	100380
	Osseospeed TX Profile	4.5/5.0, 5.0 S	Type 6	100378
	Osseospeed TX Yellow	3.0 S	Type 50	100481
Osseospeed Yellow	3.0 S	Type 50	100481	
AVINENT IMPLANT SYSTEM	Coral External	3.5	Type 30	100436
		4.1	Type 1	100353
		5.1	Type 3	100357
	Ocean External	3.5	Type 30	100436
		4.1	Type 1	100353
	Internal Hexagon Coral	3.5, 4.1	Type 32	100440
	Internal Hexagon Ocean	3.5, 4.1	Type 32	100440
	AXIS	Ceramic	4.0	Type 38
BEGO IMPLANT SYSTEMS	Semados	S3.25 - S5.5, RI 3.75 - RI 5.5	Type 26	100425
	Semados Mini-Implant	2.7, 2.9, 3.1	Type 2	100355
BIODENTA	Tissue Level	NP	Type 21	100402
	Tissue Level	RP, WP	Type 47	100478
BIOHORIZONS	Internal, Tapered Internal, Single-stage	Yellow platform	Type 32	100440
		Green platform	Type 27	100431
		Blue platform	Type 33	100442
	External	3.5	Type 30	100436

IMPLANT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
		4.0-6.0	Type 1	100353
	Laser-Lok	3.0	Type 26	100425
BIOLOK	External hex	3.45, 6.5	Type 1	100353
BIOTEC SRL	KR	3.05 platform	Type 38	100455
	KW	3.65 platform	Type 32	100440
	IM		Type 27	100431
	EN		Type 30	100436
	ER		Type 1	100353
	IR, TWO		Type 32	100440
	SR		Type 6	100378
BRAINBASE	MYTIS Arrow S series	B, BW, E, EW-Type	Type 4	100350
		C-Type	Type 47	100478
BREIDENT	blueSKY		Type 49	100480
BTI IMPLANTS	Externa®	Universal, Universal Plus	Type 1	100353
	Interna®	Universal, Universal Plus, Wide	Type 27	100431
	Tiny		Type 9	100368
BTLOCK	CV1, CV2, CV3	3.75-6.5	Type 49	100480
C-TECH	EL	3.5, 4.3, 5.1	Type 26	100425
CAMLOG	Camlog Cylinder-Line / Screw-Line	3.3	Type 28	100432
	Camlog Cylinder-Line / Screw-Line	3.8, 4.3	Type 23	100411
	Camlog Cylinder-Line / Screw-Line	5.0, 6.0	Type 24	100412
	Conelog Screw-Line	3.3	Type 65	100557
	Conelog Screw-Line	3.8, 4.3	Type 58	100541
	Conelog Screw-Line	5.0	Type 7	100380
	iSY		Type 57	100497
	Root-line Promote	3.8, 4.3	Type 23	100411
	Root-line Promote	5.0, 6.0	Type 24	100412
CHAORUM	External	3.3	Type 1	100353
	Internal	4.0, 4.8	Type 47	100478
COCHLEAR BAHA	BI 300		Type 30	100436
	BIA 400	6.0, 8.0, 10.0, 12.0	Type 55	100487
	BA 400		Type 55	100487
	Flange Fixture Hexagonal connection		Type 9	100368
COMPRESS		4.0	Type 2	100355
CORTEX	Classix, Dynamix, Saturn	3.3 - 6.0	Type 12	100359
COWELLMEDI	Atlas External	3.3	Type 2	100355
		4.1, 5.0	Type 1	100353
	Atlas Submerged	4.0	Type 21	100402
DENTAL TECH		WP	Type 47	100478
DENTAURUM IMPLANTS	tioLogic	S, M, L	Type 45	100472
DENTEGRIS	Straight & Tapered	3.3	Type 26	100425
	Straight & Tapered	3.75, 4.5, 4.75	Type 32	100440
	Straight & Tapered	5.5	Type 47	100478
DENTIS DENTAL IMPLANT SOLUTION	Internal		Type 4	100350
	Save Submerged		Type 7	100380
DENTIUM	Genoss	3.0	Type 63	100554
	Implantium	3.4, 3.8, 4.3	Type 6	100378
	Super Line	4.0 - 7.0	Type 7	100380
	Simple Line	3.8, 4.3, 4.8	Type 7	100380
	Genoss	3.6, 4.3, 5.0	Type 21	100402
DIO	External FSN	3.3	Type 2	100355
		3.75	Type 1	100353
	External FTN	5.0	Type 3	100357
	SM Submerged	3.8	Type 49	100480
		4.5, 5.3	Type 47	100478
DTI		3.5	Type 57	100497
		4.2	Type 7	100380
DYNA	Dyna Helix DC	3.2, 3.6, 4.2, 5.0	Type 27	100431
	Dyna Helix ST	3.6, 4.2, 5.0	Type 27	100431
ECKERMANN LABORATORIUM	All spiral		Type 17	100390
	External-Hex		Type 1	100353
	Triplo		Type 49	100480
EMFILS	Colosso Internal Hexagon	3.3 - 5.0	Type 11	100372
ENTIFIC	Flange fixture		Type 1	100353

IMPLANT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
EUROTEKNIKA		3.75 - 6.0	Type 1	100353
FAIR IMPLANT	FairTwo	3.5, 4.2	Type 53	100484
		5.0, 6.0	Type 6	100378
FRIALIT-2		3.4, 6.5	Type 21	100402
GC IMPLANT	GENESiO		Type 21	100402
	SETiO		Type 1	100353
	External plus	3.8 - 5.0	Type 1	100353
	External plus, Internal plus	3.0, 3.4	Type 21	100402
	Aadva Narrow		Type 38	100455
	Aadva Regular, Wide		Type 58	100541
GHIMAS	PFON EVO SERIES	3.6, 4.0	Type 38	100455
		5.0	Type 7	100380
HIOSSEN	HGII, HGIII	Mini	Type 5	100376
	HGII	Regular, Wide, Ultra Wide	Type 6	100378
	ETIII	4.0 - 5.0	Type 47	100478
	ETIII	Mini 3, 5	Type 5	100376
	ETIII	Ultra Wide 6, 7	Type 6	100378
	HSII	Regular, Wide	Type 4	100350
IDI EVOLUTION	Vite	3.75 - 4.0	Type 18	100398
IMPLANT DIRECT	ScrewPlant	3.7	Type 26	100425
	ScrewPlant	4.7, 5.7	Type 27	100431
IMTEC	Hexed-head	3.75 - 4.75	Type 1	100353
INTRA-LOCK	Internal	SQ	Type 6	100378
		Standard, Wide	Type 47	100478
ITERUM	Conical	3.75	Type 26	100425
JDENTALCARE	JDEvolution	3.7, 4.3, 5.0	Type 32	100440
		3.2	Type 38	100455
JMM	POI 2Piece	3.2, 3.7, 4.2	Type 18	100398
	POIEX HACEX	3.7	Type 15	100386
	POIEX HACEX	4.2	Type 49	100480
	POIEX HACEX	4.7, 5.2	Type 47	100478
KENTEC	SB1	3.3	Type 6	100378
	SB2	4.2	Type 6	100378
KEYSTONE DENTAL	Genesis	3.8	Type 27	100431
		4.5, 5.5, 6.5	Type 33	100442
	Renova	SDI	Type 26	100425
		RDI	Type 27	100431
	Restore	3.3	Type 30	100436
		3.75, 4.0	Type 1	100353
		5.0, 6.0	Type 3	100357
	Prima Connex	SD	Type 27	100431
		RD, WD	Type 33	100442
KLOCKNER IMPLANT SYSTEM	Essential Cone	3.5 - 5.0	Type 49	100480
	NK2, S3M, S4, S6, SK, SK2		Type 30	100436
	KL	3.7 - 4.2	Type 1	100353
LASAK	BioniQ Q-lock	Regular	Type 57	100497
	Impladent	2.9	Type 2	100355
		3.7	Type 36	100447
		5.0	Type 37	100448
LEADER	Tixos IH	3.75 - 5.5	Type 26	100425
	Tixos ITMC	3.75 - 5.5	Type 57	100497
	Implus IH	3.3	Type 26	100425
	Implus EH	3.3 - 5.0	Type 1	100353
MEDENTIKA	M-Implant	3.5-5.0	Type 38	100455
MEDICAL INSTINCT	Bone Trust	3.4 - 5.0	Type 44	100466
	Bone Trust Plus Cone	3.4, 4.0	Type 44	100466
MEDIGMA TECHNOLOGIES	Fix-A-Dent, FixTite, FixTite-S		Type 49	100480
	Fix 1		Type 17	100390
MEGAGEN	AnyOne	3.5 - 4.0	Type 7	100380
	Rescue External	7.0	Type 3	100357
	EZ Plus Internal	4.0, 4.5, 5.0	Type 7	100380
	EZ Plus Internal	Small	Type 38	100455
	EZ Plus Internal	Re/Wi	Type 7	100380
	AnyRidge		Type 27	100431

IMPLANT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
MICRODENT	Duraplant	2.2	Type 10	100370
	Microdent system		Type 1	100353
	MS, MST	3.5	Type 30	100436
	MK	5.6	Type 3	100357
	Trylogic	3.8, 4.3	Type 27	100431
	Universal	4.1, 5.1	Type 1	100353
MIS	Internal Hexagon		Type 32	100440
	Internal Octagon		Type 4	100350
	C1 (Conical Connection)	Standard 3.7	Type 49	100480
	Seven narrow	3.3	Type 63	100554
MOZO-GRAU	MG Osseous	Mini, Standard	Type 1	100353
		Maxi	Type 3	100357
	MG Inhex	Standard	Type 38	100455
		Maxi	Type 7	100380
NEMRIS	Aesthura Classic	3.75, 4.75	Type 35	100446
NEOBIOTECH	Neoplant External	4.0, 5.0	Type 1	100353
NEODENT	Alvim CM	3.5-5.0	Type 16	100388
	Alvim Ti	4.3	Type 1	100353
	Alvim II Plus	3.75-5.0	Type 21	100402
	Titamax CM	3.5-5.0	Type 16	100388
	Titamax GT	3.7-4.5	Type 9	100368
	Titamax Ti	3.75, 4.0	Type 1	100353
	Titamax Ti	4.5, 5.0	Type 3	100357
	Titamax II Plus	3.75-5.0	Type 21	100402
NEOSS	Neoss Implant	3.5-5.5	Type 18	100398
		3.25	Type 38	100455
NOBEL BIOCARE	NobelActive	NP	Type 60	100543
		NP 3.0	Type 62	100545
		RP	Type 61	100544
	Brånemark System Novum		Type 1	100353
	Brånemark System, NobelSpeedy Groovy	NP	Type 2	100355
		RP	Type 1	100353
		WP	Type 3	100357
	NobelDirect	NP	Type 8	100366
		RP	Type 9	100368
		WP	Type 9	100368
	NobelReplace, Replace Select, NobelSpeedy Replace	NP	Type 12	100359
		RP	Type 13	100361
		WP	Type 13	100361
		6.0	Type 40	100462
	Brånemark System Shorty, NobelSpeedy Shorty	RP	Type 1	100353
		WP	Type 3	100357
	NobelReplace Conical Connection	NP	Type 60	100543
		RP	Type 61	100544
	NobelReplace Platform Shift	NP 4.3	Type 12	100359
		RP 5.0, WP 6.0	Type 13	100361
NUCLEOSS	T3	3.4, 3.8, 4.2, 5.0	Type 5	100376
	T4	3.4, 3.8, 4.2	Type 5	100376
	T5	5.0	Type 7	100380
OCO BIOMEDICAL	TSI	3.25	Type 35	100446
	TSI, ERI	4.0, 5.0	Type 7	100380
	Enagage	3.2	Type 38	100455
	Enagage	3.7/3.5, 4.7/4.5	Type 26	100425
	Enagage	5.7/4.5	Type 27	100431
	Macro	6.0	Type 47	100478
OSSEOLINK		3.5-5.0	Type 21	100402
OSSTEM IMPLANTS	USII, USIII	Mini	Type 2	100355
		Regular	Type 1	100353
		Wide, Ultra Wide	Type 3	100357
	USIV	Conical Mini	Type 8	100366
	GSII	Mini	Type 5	100376
		Regular, Wide, Ultra Wide	Type 6	100378
	GSIII	Mini	Type 5	100376
		Regular	Type 6	100378
	SSI, SSII, SSIII	Regular, Wide	Type 4	100350
	TSIII SA	4.0 - 5.0	Type 47	100478
	TS III SA Fixture – NoMount	4.5	Type 47	100478

IMPLANT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
OSTEOPANT	Two piece		Type 10	100370
OSTEO-TI		3.25 - 5.5	Type 26	100425
OTICON MEDICAL	Ponto Implant		Type 9	100368
	Ponto implant Wide		Type 9	100368
PALTOP	NP		Type 38	100455
	SP		Type 49	100480
P-I BRÄNEMARK	Internal Grip Amplified	4.3, 5.1	Type 51	100482
PHILOSOPHY		3.5	Type 63	100554
	Internal Grip External Hex	3.5	Type 38	100455
	Internatl Grip External Hex	3.75-4.1	Type 18	100398
	Morse Taper	3.3, 3.75, 4.8	Type 21	100402
PHIBO	Aurea NP		Type 63	100554
	Aurea RP		Type 27	100431
	Aurea WP		Type 64	100555
	BNT Serie 3, 4, 5		Type 1	100353
	TSA Serie 3		Type 38	100455
	TSA Serie 4, 5		Type 27	100431
	TSA Advance S3		Type 38	100455
	TSA Advance S4, S5		Type 27	100431
	TSH S3, S4, S5		Type 1	100353
PLAN 1 HEALTH	Hexcel IHC, IHT, Internal Hex		Type 26	100425
	Hexcel EHC, External Hex		Type 1	100353
	GTB 3.3, 3.6, 4.3		Type 57	100497
PLATON	Type I, Type II		Type 4	100350
	Type IV Pro		Type 47	100478
	EL implant Regular 5,1		Type 47	100478
	EL implant Narrow 3,7		Type 60	100543
PRODENT ITALIA	Perform SM	3.6, 4.0	Type 49	100480
		4.5, 5.0	Type 7	100380
REVOIS		3.8, 4.3, 5.0	Type 26	100425
SDS	SDS 2.0	3.8	Type 69	100589
SIC IMPLANT SYSTEM	SIC ace	3.4, 4.0, 5.0	Type 44	100466
S.I.N. IMPLANT SYSTEM	External Hex	3.4, 4.1, 5.0	Type 1	100353
	Internal Hex	3.8, 4.5, 5.5	Type 15	100386
	Strong SW HI	4.5	Type 48	100479
	Strong SW Morse	4.5	Type 53	100484
SERVO-DENTAL	Trias	3.8	Type 6	100378
SIC IMPLANT SYSTEM	SIC ace	3.4, 4.0, 5.0	Type 44	100466
SOUTHERN IMPLANTS	External Hex		Type 1	100353
	Extra-Oral		Type 1	100353
	IT-Connection		Type 4	100350
	MAX 7.0, 8.0, 9.0		Type 1	100353
	MAX-IT 7.0, 8.0, 9.0		Type 7	100380
	Tri-MAX 7.0		Type 24	100412
	Tri-MAX 8.0, 9.0		Type 40	100462
	Tri-nex 3.5		Type 12	100359
	Tri-nex 4.3, 5.0		Type 13	100361
STERNGOLD DENTAL	Stern IC Internal Octagon	Regular Neck 4.1	Type 4	100350
	External Hexagon	3.75, 4.0, 5.0 RP	Type 1	100353
STRAUMANN DENTAL	Bone Level Narrow CrossFit	3.3	Type 53	100484
IMPLANT SYSTEM	Bone Level Regular CrossFit	4.1, 4.8	Type 54	100485
	Extra Oral		Type 1	100353
	Ortho/Palatal		Type 9	100368
	Roxolid: Standard, Standard Plus, TE	Regular Neck 3.3	Type 4	100350
	Roxolid: Bone Level Narrow CrossFit	3.3	Type 53	100484
	Roxolid: Bone Level Regular CrossFit	4.1, 4.8	Type 54	100485
	Roxolid: Standard Standard Plus, Narrow Neck		Type 17	100390
	Roxolid: Standard, StandardPlus, Narrow Neck/Wide Neck		Type 4	100350
	Standard, Standard Plus, TE	Narrow Neck	Type 17	100390
		Regular Neck/Wide Neck	Type 4	100350
	NNC	3.3	Type 68	100560
SWEDEN & MARTINA	Kohno - new platform	3.8	Type 32	100440
	Kohno - old platform	3.8	Type 21	100402
	Kohno - old + new platform	4.25, 5.0	Type 7	100380
	Premium - old platform	3.3, 3.8	Type 21	100402

IMPLANT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
	Premium - old + new platform	4.25, 5.0	Type 7	100380
	Global		Type 48	100479
SYBRON	Endopore External	4.1, 5.0	Type 1	100353
	Endopore External	3.5	Type 30	100436
	Entegra External	3.25, 4.0	Type 1	100353
	Pitt-Easy	4.0, 4.9	Type 24	100412
	Sybron PRO XRT	3.3, 3.5, 4.1	Type 6	100378
TA-DENT		3.3	Type 38	100455
		4.8	Type 6	100378
		5.6	Type 49	100480
TBR	External Hexagon		Type 1	100353
	Internal Octagon		Type 18	100398
	M-implant, Z1-M implant	3.2 / 3.9 / 4.7	Type 21	100402
TFI SYSTEMS		3.3 - 5.0	Type 26	100425
TEKKA	In-Kone, T-Quest		Type 49	100480
	Krestal, Progress		Type 6	100378
	TwinKon		Type 17	100390
THOMMEN MEDICAL	SPI System	platform 3.5	Type 22	100404
		platform 4.0	Type 21	100402
		platform 4.5 - 6.0	Type 15	100386
TIDAL SPIRAL IMPLANT SYST.		3.75, 4.2, 5.0, 6.0	Type 27	100431
TRI DENTAL IMPLANTS	Tri-Vent	3.5 platform	Type 32	100440
U-IMPL	Internal Hexagon		Type 32	100440
	External Hexagon		Type 1	100353
WARANTEC	Oneplant internal tapered		Type 6	100378
	Oneplant implant	4,3	Type 24	100412
WIELAND DENTAL IMPLANTS		3.5, 4.3, 5.0	Type 45	100472
XIVE		3.0	Type 22	100404
		3.4	Type 21	100402
		3.8 - 5-5	Type 45	100472
		TG	Type 8	100366
ZIMMER DENTAL	AdVent	4.5	Type 35	100446
		5.7	Type 34	100444
	Screw-Vent (also Tapered)	3.5 (3.5, 4.1)	Type 32	100440
		4.5 (4.7)	Type 27	100431
		5.7 (6.0)	Type 33	100442
	Spline	3.25 - 5.0	Type 30	100436
	SwissPlus (also Tapered)	platform 3.8	Type 31	100438
		platform 4.8	Type 4	100350

ABUTMENT SYSTEM	MODEL	SIZE	SMARTPEG	ARTICLE NO.
COCHLEAR BAHA	BA 210		Type 55	100487
	BA 300		Type 55	100487
NEOSS	Access Abutment		Type 52	100483
NOBEL BIOCARE	Brånemark System MUA	NP, RP	Type A3	100363
OTICON MEDICAL	Ponto Implant Abutment	4.0	Type 55	100487
THOMMEN MEDICAL	SPI SYSTEM Bar Abutment		Type 36	100447
	VARIOMulti Abutment		Type 25	100415

Please note that this document is updated continuously!
The latest version can always be found at www.osstell.com

Snap this QR code to view the latest SmartPeg reference list on your Smartphone.



25016-56