




BEGO ON TISSUE LEVEL



BEGO SEMADOS® S-IMPLANTATE TG
 QUICK & FUNCTIONAL · UNIVERSAL & SUCCESSFUL
 RELIABLE & PROVEN · VARIED & SOPHISTICATED
 100% GERMAN QUALITY PRODUCT
 SIMPLE & INTUITIVE · QUICK & FUNCTIONAL
 UNIVERSAL & SUCCESSFUL · RELIABLE & PROVEN
 VARIED & SOPHISTICATED · SIMPLE & INTUITIVE

Reaching new heights

Tailored to your patient down to the finest detail. Various implant systems, supplemented by individual CAD/CAM fabricated prosthetic components and surgical guides as well as the BEGO Multi^{Plus} concept—a system that provides patients with a perfectly fitting prosthesis in just one day. The exceptionally high technical standards ensure a high degree of satisfaction—for both dental implantologists and patients.

Your benefits

- Implant interface at gingival height
- High machined implant shoulder
- One-stage implantation technique is possible
- Same implant-abutment interface as for BEGO Semados® SC/SCX/RS/R SX/RI implant systems
- Stable connection with no micro-movements during functional masticatory loading (up to 150 N)
- Use in combination with BEGO Semados® S-Line Tray^{Plus}
- For the most stringent quality requirements—100% German development and manufacturing

Technical details

- Self-tapping implant with Platform Switch and 3 mm machined implant shoulder
- Manufactured from grade 4 pure titanium with ultrapure homogeneous TiPure^{Plus} surface
- Diameter: 3,75 · 4,1 mm
- Lengths: 10 · 13 mm (osseous nominal size + 2,5 mm gingival area)
 - S-Implant TG L10, total length = 12,5 mm
 - S-Implantate TG L13, total length = 15,5 mm
- Gap-free fit of the abutments thanks to the internal taper

PRODUCT OVERVIEW

Distributor contact data

BEGO Semados® S-Implants TG (PS)

Art. no.	Color code	Product	Quantity
BEGO Semados® S-Implants TG 3,75			
58330	■	S 3,75 TG L 10	—
58331	■	S 3,75 TG L 13	—
BEGO Semados® Mini-Implants 4,1			
58333	■	S 4,1 TG L 10	—
58334	■	S 4,1 TG L 13	—